allow you to slide one finger between the strap and your leg. If the brace slips during use and you need additional stabilisation, use the additional strap number 5. Before putting the brace on, fasten the additional strap around the popliteal fossa (Fig. 5), so that the loops on the belt are on the back side of the limb. Proceed to putting on the orthosis, thread strap number 1 through the belt loops number 5 (Fig. 6). Fasten the straps as shown in Fig. 4.

Size chart (left and right foot braces are available):

Size	S	M	L	XL
Thigh circumference [cm]	39-47	47-53	53-59	59-67
Knee circumference (cm)	33-35	35-38	38-43	43-48
Calf circumference (cm)	30-35	35-40	40-47	47-50

Thigh circumference measurement – 15 cm above the centre of the kneecap

Knee circumference measurement - at the centre of the kneecap

Calf circumference measurement - 15 cm below the centre of the kneecap

Cleaning instructions:

Clean the frame components with a damp cloth. Wipe the fabric padding elements with a damp cloth or hand wash at low temperatures (up to 30°) with very mild detergents, without softeners. Do not chlorinate, do not dry-clean, do not wring, do not spin, squeeze gently before drying. Dry away from artificial heat sources. Dry unfolded. Do not iron.

Composition:

PA, AL, PE, PU, PP.

Intended user:


The medical device is intended for regular (non-professional) users, as well as for professional users.

WARNINGS:

1. The device should not be used by children without adult supervision.
2. Do not use the device on uncleaned skin.
3. Do not use the device if it has been damaged or deformed.
4. Do not attempt to repair the device yourself.
5. Do not use the device made for another patient or use it in a manner for which it is not intended.
6. Do not use the device in case of prolonged immobilisation of the limb in one position.
7. The product is intended to be worn directly on the leg.
8. Do not use the brace on skin covered with warming gels or ointments.
9. Failure to maintain proper personal hygiene or the use of inappropriate detergents may cause abrasions, increased sweating, or the development of bacterial flora.
10. In case of any problems related to the use of the product, contact your doctor.
11. Any serious incident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IT IS RECOMMENDED TO CONSULT
A SPECIALIST BEFORE USING THE DEVICE

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EN

JOINTCRAFT 3D KNEE BRACE
WITH ADJUSTABLE RANGE OF MOTION

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Device description:

The knee brace has a four-axis mechanism with adjustable angular range of flexion (10°, 20°, 30°, 45°, 60°, 90°) and extension (10°, 20°, 30°) as well as immobilisation at 0°. It has soft upholstery which, thanks to its non-slip properties, prevents the brace from slipping off the limb. Stabilising straps ensure the correct positioning of the brace. For an even better adjustment of the brace to the patient's limb, replaceable condyle pads (under the mechanism) of various thicknesses are included, which prevent skin irritation and chafing. The brace comes in two versions - personalised and universal. In the personalised version, three-dimensional scanning technology of the patient's limb is used, which enables an accurate reflection of the shape. It provides the starting model for the development of a 3D model of the device, which is then printed from materials designed for this purpose.

Indications:

The brace is intended for people with knee joint instability of various aetiologies: it can be used for ACL, PCL, MCL, LCL ligament injuries, and as a form of conservative treatment among patients with osteoarthritis of the knee, and it helps in the rehabilitation process after surgery as well as in the return to full fitness.

Contraindications:

Varicose veins of the lower limb, acute or chronic inflammation of the knee joint, severe deformations, post-traumatic or post-operative conditions requiring complete immobilisation.

Target patient groups:

Patients with knee joint instability of various aetiologies.

Disposal of the medical device:

The device must be disposed of following the applicable local environmental regulations.

Choosing your brace:

1. The size of a standard brace should be selected based on the dimensions shown in the size chart.
2. A personalised brace is a product made to order based on the individual measurements of the patient.
3. If the device becomes too loose or changes its position, it must be refitted starting from the first step in the section „How to put on the brace”.

How to put the orthosis on:

1. It is recommended that the brace is fitted for the first time by qualified medical personnel.
2. Before fitting the brace, set the correct range of motion restrictions as recommended by your doctor or physiotherapist. Make sure the settings on both dials are symmetrical.
3. All restrictors are marked with a letter indicating the side they should be used on, i.e. L-left or R-right, the angle to which the movement is restricted and which movement is restricted, i.e. F (flexion) and E (extension) (Fig. 3).
4. Place the chosen extension angle restrictor on the front of the dial and screw it in place with the screw provided (Fig. 1).
5. Repeat on the back of the dial, placing a restrictor with the mark of the chosen flexion angle range (Fig. 2).
6. Undo all straps.
7. Fit the brace so that the frame components rest on the front of the thigh and on the calf, and the side splints are located symmetrically on both sides of the knee.
8. The axis of the dial should be at the level of the upper edge of the patella (Fig. 4). Check that the dial fits snugly on both sides of the knee joint. If not, use the appropriate condyle pads to adjust the dial.
9. Fasten the straps in a strict sequence (Fig. 4): 1. upper tibial strap. 2. lower tibial strap. 3. upper femoral strap. 4. lower femoral strap.
10. The stabilising straps should not be over-tightened. The strap should be loose enough to